

Product Datasheet

PTAFR antibody (FITC) (orb25463)

Description	PTAFR antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Platelet-activating factor receptor protein (297-342AA)
Target	PTAFR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25105
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.