



Product Datasheet

RANBP2 antibody (orb25511)

Description	RANBP2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human E3 SUMO-protein ligase RanBP2 protein (2601-2802AA)
Target	RANBP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49792
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue showing RANBP2 expression.



Immunohistochemical staining of human testis tissue showing RANBP2 expression.