

Product Datasheet

RBM6 antibody (orb25536)

Description	RBM6 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human RNA-binding protein 6 protein (300-535AA)
Target	RBM6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P78332
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.