


Product Datasheet

TRAPPC10 antibody (orb255457)

Description	TRAPPC10 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P48553
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of U2OS
cell...