


Product Datasheet

STAB1 Antibody (orb255657)

Description	STAB1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9NY15
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Co...



Immunohistochemical staining of Human Th...