



Product Datasheet

DDK1 Antibody (orb255829)

Description	DDK1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O94907
Expiration Date	12 months from date of receipt.



Immunohistochemical
analysis of Human
li...



Immunohistochemical
analysis of Human
br...