



Product Datasheet

Apolipoprotein CIII antibody (orb256341)

Description	Apolipoprotein CIII antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	Human Human Apolipoprotein CIII recombinant protein
Target	Apolipoprotein CIII
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	After reconstitution store at 4°C.for up to two weeks. For long term storage, aliquot and store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	9 kDa
Uniprot ID	P02656, P33622, P06759
NCBI	000040, 000031
Entrez	345
Dilution Range	IHC-P: 1:200
Expiration Date	12 months from date of receipt.



IHC-P image of rat skin tissue using ant...



Immunohistochemical staining of paraffin...