



Product Datasheet

RANTES antibody (orb256342)


Description	RANTES antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	Human recombinant RANTES protein. Please contact us for the exact immunogen sequence. The protein is available as orb334647
Target	RANTES
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	After reconstitution store at 4°C.for up to two weeks. For long term storage, aliquot and store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	10 kDa
Uniprot ID	P13501
Dilution Range	IF/ICC: 1:200, IHC-P: 1:500, WB: 1:200-1000
Expiration Date	12 months from date of receipt.



IF analysis of human
hela cells tissue u...



IHC-P image of pig
small intestine tissu...



Immunohistochemical
staining of pig larg...