

## Product Datasheet

### Leptin Receptor antibody (orb256343)

## Description

Leptin Receptor antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

### Immunogen

KLH conjugated synthetic peptide derived from human Leptin Receptor. Please contact us for the exact immunogen sequence. The peptide is available as orb374876.

### Target

Leptin Receptor

### Form/Appearance

10 mM PBS, 0.02% sodium azide

### Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

### Storage

After reconstitution store at 4°C.for up to two weeks. For long term storage, aliquot and store at -20°C. Avoid freeze/thaw cycles.

### Note

For research use only

### Isotype

IgG

### Clonality

Polyclonal

### Purity

Polyclonal antibodies are purified by peptide affinity chromatography

### MW

132 kDa

### Uniprot ID

[Q6PDU4](#), [Q9CQ74](#), [O95214](#)

### NCBI

[001121680](#), [056159](#)

### Entrez


[23484](#)

### Dilution Range

WB: 1:200-600

### Expiration Date

12 months from date of receipt.



WB analysis of human breast cancer 24 (l...