

Product Datasheet

MCP3 antibody (orb256344)

Description

MCP3 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P

Applications

Immunogen

KLH conjugated synthetic peptide derived from human MCP3. Please contact us for the exact immunogen sequence. The peptide is available as orb374877.

Target

MCP3

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

After reconstitution store at 4°C.for up to two weeks. For long term storage, aliquot and store at -20°C. Avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

11 kDa

Uniprot ID

[P80098](#), [Q9QXY8](#), [Q03366](#)

NCBI

[038682](#), [013654](#), [31](#), [11](#)

Entrez


[20306](#)

Dilution Range


IF/ICC: 1:200, IHC-P: 1:200-500

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunofluorescence analysis of mouse ski...



IHC-P staining of pig kidney tissue usin...