

Product Datasheet

MIP3 beta antibody (orb256345)

Description

MIP3 beta antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P

Applications

Immunogen

KLH conjugated synthetic peptide derived from human MIP3 beta. Please contact us for the exact immunogen sequence. The peptide is available as orb374878.

Target

MIP3 beta

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

After reconstitution store at 4°C.for up to two weeks. For long term storage, aliquot and store at -20°C. Avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

9.2 kDa

Uniprot ID

[Q99731](#), [O70460](#), [D3Z184](#)

NCBI

[005064](#), [005055](#)

Entrez



[6363](#)

Dilution Range


IF/ICC: 1:200, IHC-P: 1:200

Expiration Date


12 months from date of receipt.



IHC-P staining of mouse lymph node tissue...



Immunohistochemical staining of rat colo...



IF analysis of mouse skin tissue using a...