



Product Datasheet

Transthyretin antibody (orb256346)


Description	Transthyretin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant human Transthyretin protein. The protein is available as (orb222833) was used as immunogen.
Target	Transthyretin
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	After reconstitution store at 4°C. for up to two weeks. For long term storage, aliquot and store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	16 kDa
Uniprot ID	P02766, P07309
NCBI	000362, 000371
Entrez	7276
Dilution Range	WB: 1:200-1000
Expiration Date	12 months from date of receipt.



WB analysis of rat liver (lane 1), rat t...



WB analysis of Transthyretin-his recombi...



Western blot analysis of rat liver (lane...)