



Product Datasheet

c-myc epitope tag antibody (orb256348)

Description	Mouse monoclonal antibody to c-myc tag
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P, IP, WB
Immunogen	A synthetic peptide corresponding to aa 408-439 from C-terminus of human c-myc.
Target	c-myc epitope tag
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P01106
Expiration Date	12 months from date of receipt.



ELISA of
anti-epitope
tag
antibodies
on ...



WB analysis
of of
HEK293
spiked (B)
a...