

Product Datasheet

c-myc tag antibody (orb256351)

Description	Human monoclonal antibody to c-myc tag
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P, IP, WB
Immunogen	A synthetic peptide corresponding to aa 408-439 from C-terminus of human c-myc.
Target	c-myc epitope tag
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	Fab fragment
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P01106
Expiration Date	12 months from date of receipt.