

Product Datasheet

Fluorescein antibody (orb256391)

Description	Mouse monoclonal antibody to Fluorescein
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	Control
Immunogen	Fluorescein I-bovine serum albumin (FI-BSA).
Target	Fluorescein
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody is not conjugated to fluorescein, but was raised against fluorescein. This antibody can be used as an isotype control in both in vitro and in vivo applications.
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.