

Product Datasheet

hapten 4-hydroxy-3-NP acetyl antibody (orb256409)

Description	Mouse monoclonal antibody to hapten 4-hydroxy-3-NP acetyl
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	Control
Immunogen	hapten-protein conjugated NP12-chicken gammaglobulin (CG).
Target	hapten 4-hydroxy-3-nitrophenyl acetyl (NP)
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2c
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.