

Product Datasheet

Beta galactosidase antibody (orb256416)

Description	Mouse monoclonal antibody to Beta galactosidase
Species/Host	Human
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Beta-galactosidase.
Target	Beta-galactosidase
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	Fab fragment
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P00722
Expiration Date	12 months from date of receipt.