


Product Datasheet

V5 tag antibody (orb256432)

Description	Mouse monoclonal antibody to V5 tag
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IP, RIA, WB
Immunogen	Pk (V5) Epitope Tag (GKPIPPLLGLDST).
Target	V5 epitope tag
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.



ELISA of anti-epitope tag antibodies on

...