

## Product Datasheet

### V5 tag antibody (orb256441)

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit monoclonal antibody to V5 tag                                             |
| <b>Species/Host</b>        | Mouse                                                                            |
| <b>Reactivity</b>          | Virus                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                     |
| <b>Tested Applications</b> | ELISA, IF, IHC-Fr, IP, RIA, WB                                                   |
| <b>Immunogen</b>           | Pk (V5) Epitope Tag (GKPIPPLLGLDST).                                             |
| <b>Target</b>              | V5 epitope tag                                                                   |
| <b>Preservatives</b>       | PBS with 0.02% Proclin 300.                                                      |
| <b>Concentration</b>       | 1 mg/ml                                                                          |
| <b>Storage</b>             | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. |
| <b>Note</b>                | For research use only                                                            |
| <b>Isotype</b>             | IgG                                                                              |
| <b>Clonality</b>           | Monoclonal                                                                       |
| <b>Purity</b>              | Purified                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                  |