



Product Datasheet

HER2 antibody (orb256519)

Description	Rabbit polyclonal antibody to HER2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y1221/Y1222 of human HER2 protein. The exact sequence is proprietary.
Target	ERBB2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P04626
Entrez	2064
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IP: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of HER2 (Phospho-Y1221)



Immunohistochemical analysis of HER2 (Phospho-Y1221)