



Product Datasheet

HER3 (phospho-Y1328) antibody (orb256520)


Description	Rabbit polyclonal antibody to HER3 (phospho-Y1328)
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y1328 of human HER3 protein. The exact sequence is proprietary.
Target	ERBB3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P21860, Q61526, Q62799
Entrez	2065, 29496, 13867
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of HER3 (Phospho-Y1328) antibody.



Immunohistochemical analysis of HER3 (Phospho-Y1328) antibody.



Immunofluorescent analysis of HER3 (Phospho-Y1328) antibody.