

Product Datasheet

Tissue Factor (phospho-S290) antibody (orb256529)

Description

Rabbit polyclonal antibody to Tissue Factor

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

IH, WB

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S290 of human CD142 protein. The exact sequence is proprietary.

Target

F3

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[P13726](#)

Entrez


[2152](#)

Dilution Range

WB: 1:500-1000

Expiration Date

12 months from date of receipt.



Western blot analysis of CD142 (Phospho-...



Immunohistochemical analysis of CD142 (P...