



Product Datasheet

FRS2 (phospho-Y436) antibody (orb256541)

Description	Rabbit polyclonal antibody to FRS2 (phospho-Y436)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y436 of human FRS2 protein. The exact sequence is proprietary.
Target	FRS2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8C180 , Q8WU20
Entrez	10818 , 327826
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Western blot analysis of FRS2 (Phospho-Y436) antibody



Immunohistochemical analysis of FRS2 (Phospho-Y436) antibody