


Product Datasheet

Frizzled 1/2/7 antibody (orb256544)

Description	Rabbit polyclonal antibody to Frizzled 1/2/7
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human Frizzled 1/2/7. The exact sequence is proprietary.
Target	FZD1; FZD2; FZD7
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q61090 , Q9UP38 , O75084 , Q14332
Entrez	14369 , 8324
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Frizzled 1/2/7 ...



Immunofluorescent analysis of Frizzled 1...