

Product Datasheet

GATA1 antibody (orb256551)

Description	Rabbit polyclonal antibody to GATA1
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Primate, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ChIP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human GATA1. The exact sequence is proprietary.
Target	GATA1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P17679 , P43429 , P15976
Entrez	100911167 , 25172 , 2623 , 14460
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of GATA1 expression...