



Product Datasheet

HAND1 antibody (orb256583)

Description	Rabbit polyclonal antibody to HAND1
Species/Host	Rabbit
Reactivity	Human, Porcine, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human HAND1. The exact sequence is proprietary.
Target	HAND1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P97832 , O96004 , Q64279
Entrez	59112 , 15110 , 9421
Dilution Range	WB: 1:500-1000, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of HAND1 expression...



Immunofluorescent analysis of HAND1 staining...