

Product Datasheet

Histone Deacetylase 6 (phospho-S22) antibody (orb256587)

Description

Rabbit polyclonal antibody to Histone

Species/Host

Rabbit

Reactivity

Human, Primate

Conjugation

Unconjugated

Tested Applications

IF, IH, WB

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S22 of human Histone Deacetylase 6 protein. The exact sequence is proprietary.

Target

HDAC6

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[Q9Z2V5](#), [Q9UBN7](#)

Entrez

[15185](#), [10013](#)

Dilution Range

WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500

Expiration Date

12 months from date of receipt.

kDa A


180

130


95

72

Western blot analysis of Histone Deacety...



Immunohistochemical analysis of Histone ...



Immunofluorescent analysis of Histone De...