



Product Datasheet

hnRNP A2/B1 antibody (orb256605)


Description	Rabbit polyclonal antibody to hnRNP A2/B1
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human hnRNP A2/B1. The exact sequence is proprietary.
Target	HNRNPA2B1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	A7VJC2 , O88569 , P22626
Entrez	3181 , 362361 , 53379 , 102642938
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of hnRNP A2/B1 exp...



Immunohistochemical analysis of hnRNP A2...



Immunofluorescent analysis of hnRNP A2/B...