



Product Datasheet

MEKK1 (phospho-T1402) antibody (orb256669)

Description	Rabbit polyclonal antibody to MEKK1 (phospho-T1402)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T1402 of human MEKK1 protein. The exact sequence is proprietary.
Target	MAP3K1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q13233 , P53349 , Q62925
Entrez	116667 , 4214
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Western blot analysis of MEKK1 (Phospho-T1402) antibody.



Immunohistochemical analysis of MEKK1 (Phospho-T1402) antibody.