



Product Datasheet

Tau antibody (orb256674)


Description	Rat polyclonal antibody to Tau having MW of 82 kD that helps stabilize microtubules in our body. On ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S673 of human TAU protein. The exact sequence is proprietary.
Target	MAPT
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P10636, P19332, P10637
Entrez	17762, 4137
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of TAU (Phospho-S673)



Immunohistochemical analysis of TAU (Phospho-S673)



Immunofluorescent analysis of TAU (Phospho-S673)