



Product Datasheet

NFAT4 (phospho-S165) antibody (orb256716)


Description	Rabbit polyclonal antibody to NFAT4 (phospho-S165)
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S165 of human NFAT4 protein. The exact sequence is proprietary.
Target	NFATC3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q12968 , P97305
Entrez	4775
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of NFAT4 (Phospho-S165) antibody



Immunohistochemical analysis of NFAT4 (Phospho-S165) antibody



Immunofluorescent analysis of NFAT4 (Phospho-S165) antibody