


Product Datasheet


NF-kappaB p105/p50 (phospho-S337) antibody (orb256719)

Description


Rabbit polyclonal antibody to NF-kappaB



Western blot analysis of NF-kappaB p105/...



Immunohistochemical analysis of NF-kappa...



Immunofluorescent analysis of NF-kappaB ...

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested Applications

IF, IH, WB

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S337 of human NF-kappaB p105/p50 protein. The exact sequence is proprietary.

Target

NFKB1

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[Q63369](#), [P25799](#), [P19838](#)

Entrez

[18033](#), [4790](#)

Dilution Range

WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500

Expiration Date

12 months from date of receipt.