

Product Datasheet

E-selectin antibody (orb25673)

Description	E-selectin antibody
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Pig E-selectin protein (23-429AA)
Target	SELE
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P98110
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.