



Product Datasheet

PAK1 (phospho-T212) antibody (orb256740)


Description	Rabbit polyclonal antibody to PAK1 (phospho-T212)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T212 of human PAK1 protein. The exact sequence is proprietary.
Target	PAK1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P35465, O88643, Q13153
Entrez	29431, 5058
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of PAK1 (Phospho-T212)



Immunohistochemical analysis of PAK1 (Phospho-T212)



Immunofluorescent analysis of PAK1 (Phospho-T212)