

Product Datasheet

GABBR1 peptide (orb256836)

Description	The peptide is used to anti GABBR1 antibody reactivity.
Reactivity	Human, Primate, Rat
Tested Applications	Blocking
Target	GABBR1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256546 (GABBR1 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9WV18 , Q9UBS5 , Q9Z0U4
Entrez	54393 , 2550 , 81657
Expiration Date	12 months from date of receipt.