

Product Datasheet

Histone Deacetylase 5/9 peptide (orb256876)

Description	The peptide is used to anti Histone Deacetylase 5/9 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	HDAC5; HDAC9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256586 (Histone Deacetylase 5/9 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q99N13 , Q9UKV0 , Q9Z2V6 , Q9UQL6
Entrez	9734 , 79221 , 10014
Expiration Date	12 months from date of receipt.