

Product Datasheet

CD18 peptide (orb256927)

Description	The peptide is used to anti CD18 antibody reactivity.
Reactivity	Human, Primate
Tested Applications	Blocking
Target	ITGB2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256637 (CD18 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05107, P11835
Entrez	16414, 3689
Expiration Date	12 months from date of receipt.