

Product Datasheet

Kir6.2 peptide (orb256932)

Description	The peptide is used to anti Kir6.2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	KCNJ11
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256642 (Kir6.2 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q61743, Q14654, P70673
Entrez	16514, 3767, 83535
Expiration Date	12 months from date of receipt.