

Product Datasheet

KLF10/11 peptide (orb256936)

Description	The peptide is used to anti KLF10/11 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	KLF10; KLF11
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256646 (KLF10/11 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O89091 , Q8K1S5 , Q13118 , O14901 , O08876
Entrez	7071 , 21847 , 194655 , 81813 , 8462
Expiration Date	12 months from date of receipt.