

Product Datasheet

MYL3 peptide (orb256992)

Description	The peptide is used to anti MYL3 antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	MYL3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256702 (MYL3 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P08590
Entrez	4634
Expiration Date	12 months from date of receipt.