

Product Datasheet

p47 phox (phospho-S359) peptide (orb256998)

Description	The peptide is used to anti p47 phox (phospho-S359) antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	NCF1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256708 (p47 phox (phospho-S359) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P14598
Entrez	653361
Expiration Date	12 months from date of receipt.