

Product Datasheet

p40 phox peptide (orb256999)

Description	The peptide is used to anti p40 phox antibody reactivity.
Reactivity	Human, Porcine, Rabbit
Tested Applications	Blocking
Target	NCF4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256709 (p40 phox antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P97369, Q15080
Entrez	4689, 17972
Expiration Date	12 months from date of receipt.