

## Product Datasheet

### pro-NGF beta peptide (orb257011)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | The peptide is used to anti pro-NGF beta antibody reactivity.  |
| <b>Reactivity</b>          | Human, Mouse, Rat  |
| <b>Tested Applications</b> | Blocking   |
| <b>Target</b>              | NGF  |
| <b>Form/Appearance</b>     | Lyophilized powder   |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.   |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb256721 (pro-NGF beta antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| <b>Purity</b>              | > 85%  |
| <b>Source</b>              | Synthetic  |
| <b>Uniprot ID</b>          | <b>P01138, P01139</b>  |
| <b>Entrez</b>              | <b>18049, 4803</b>   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |