

Product Datasheet

Kappa Opioid R peptide (orb257027)

Description	The peptide is used to anti Kappa Opioid R antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	OPRK1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256737 (Kappa Opioid R antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P41145, P33534, P34975
Entrez	4986, 18387, 29335
Expiration Date	12 months from date of receipt.