

Product Datasheet

PINX1 peptide (orb257042)

Description	The peptide is used to anti PINX1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PINX1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256752 (PINX1 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9CZX5 , A4L691 , Q96BK5
Entrez	305963 , 72400 , 54984
Expiration Date	12 months from date of receipt.