



Product Datasheet

CD82 antibody (orb257047)


Description	Mouse monoclonal antibody to CD82
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, ICC, IHC-P, IP, WB
Immunogen	C91/PL (human HTLV-1+ T cell line)
Target	CD82
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P27701
Entrez	3732
Dilution Range	Functional application: In human MOLT-4 T-cell line the antibody C33 inhibits syncytium formation induced by coculture with human T-cell leukemia virus type 1 (HTLV-1)-positive human T-cell lines. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human CD82 positive lympho...



Flow cytometry surface staining pattern ...



Immunohistochemistry staining of human t...