

Product Datasheet

CD105 antibody (FITC) (orb257051)

Description

CD105 antibody (FITC)

Species/Host

Rabbit

Reactivity

Bovine, Canine, Equine, Gallus, Human, Mouse, Porcine, Rat

Conjugation

FITC

Tested

ICC, IF

Applications
Immunogen

Synthetic peptide (please contact us for details)

Target

CD105 (FITC)

Form/Appearance

Liquid: PBS, 2 mM Sodium Azide, 30% glycerol

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

70 kDa

Uniprot ID
P17813, P37176, Q63961
NCBI
NM_000118.31, NP_000109.11, NP_000109.1, NM_000118.3
Entrez
2022
Dilution Range

IF/ICC:1:50-200 (based on 0.5 mg/ml)

Expiration Date

12 months from date of receipt.