

Product Datasheet

Endothelin 1 antibody (FITC) (orb257056)

Description

Endothelin 1 antibody (FITC)

Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat, Sheep
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Synthetic peptide (please contact us for details)
Target	Endothelin 1 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	24 kDa
Uniprot ID	P05305 , P22388 , P22387
NCBI	NP_001161791.1 , NM_001168319.1
Entrez	1906
Dilution Range	IF/ICC:1:50-200 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.