

Product Datasheet

GM1(GS) antibody (FITC) (orb257057)

Description

GM1(GS) antibody (FITC)

Species/Host

Rabbit

Reactivity

Bovine, Canine, Equine, Goat, Guinea pig, Human, Mouse, Porcine, Rat

Conjugation

FITC

Tested

ICC, IF

Applications
Immunogen

Synthetic peptide (please contact us for details)

Target

GM1(GS) (FITC)

Form/Appearance

Liquid: PBS, 2 mM Sodium Azide, 30% glycerol

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

1.57 kDa

Uniprot ID
P38438, P37173, Q62312
NCBI
NM_003242.51, NP_003233.41, NM_003242.5, NP_003233.4
Entrez
7048
Dilution Range

IF/ICC:1:50-200 (based on 0.5 mg/ml)

Expiration Date

12 months from date of receipt.