

## Product Datasheet

### Integrin beta 7 (FITC) antibody (orb257058)

|                            |                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Integrin beta 7 (FITC) antibody                                                                                |
| <b>Species/Host</b>        | Rabbit                                                                                                         |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                              |
| <b>Conjugation</b>         | FITC                                                                                                           |
| <b>Tested Applications</b> | ICC, IF                                                                                                        |
| <b>Immunogen</b>           | Synthetic peptide (please contact us for details)                                                              |
| <b>Target</b>              | Integrin beta 7 (FITC)                                                                                         |
| <b>Form/Appearance</b>     | Liquid: PBS, 2 mM Sodium Azide, 30% glycerol                                                                   |
| <b>Concentration</b>       | - 100 µg (in 200 µl): 0.5 mg/ml                                                                                |
| <b>Storage</b>             | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                                          |
| <b>Isotype</b>             | IgG                                                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                                                     |
| <b>Purity</b>              | Polyclonal antibodies are purified by peptide affinity chromatography                                          |
| <b>MW</b>                  | 85 kDa                                                                                                         |
| <b>Uniprot ID</b>          | <b>P26010, Q4G029, P26011</b>                                                                                  |
| <b>NCBI</b>                | <b>NM_000889.2, NP_000880.1</b>                                                                                |
| <b>Entrez</b>              | <b>3695</b>                                                                                                    |
| <b>Dilution Range</b>      | IF/ICC:1:50-200 (based on 0.5 mg/ml)                                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                |